

Dr Emilie Marhuenda PhD

School of Engineering and Materials Science
Queen Mary University of London
Mile End Road
London E1 4NS

tel: +44 (0)20 7882 8732
email: e.marhuenda@qmul.ac.uk web: www.sems.qmul.ac.uk/e.marhuenda

2024

Sox9 Accelerates Vascular Aging by Regulating Extracellular Matrix Composition and Stiffness.

Faleeva M, Ahmad S, Theofilatos K, Lynham S, Watson G, Whitehead M, Marhuenda E, Iskratsch T, Cox S and Shanahan CM. *Circulation Research* vol. 134, (3) 307-324. Wolters Kluwer.

2023

Hypertensive Pressure Mechanosensing Alone Triggers Lipid Droplet Accumulation and Transdifferentiation of Vascular Smooth Muscle Cells to Foam Cells.

Swiatlowska P, Tipping W, Marhuenda E, Severi P, Fomin V, Yang Z, Xiao Q, Graham D, Shanahan C and Iskratsch T. *Advanced Science* vol. 11, (9). Wiley.

Recent Advances in Electrospun Fibers for Biological Applications.

Fromager B, Marhuenda E, Louis B, Bakalara N, Cambedouzou J and Cornu D. *Macromolâ€”a Journal of Macromolecular Research* vol. 3, (3) 569-613. Mdpi.

Lem2 is essential for cardiac development by maintaining nuclear integrity.

Ross JA, Arcos-Villacis N, Battey E, Boogerd C, Orellana CA, Marhuenda E, Swiatlowska P, Hodzic D, Prin F, Mohun T, Catibog N, Tapia O, Gerace L, Iskratsch T, Shah AM and Stroud MJ. *Cardiovascular Research* vol. 119, (11) 2074-2088. Oxford University Press (Oup).

2022

Photoelectrochemical imaging of single cardiomyocytes and monitoring of their action potentials through contact force manipulation of organoids.

Jacques R, Zhou B, Marhuenda E, Gorecki J, Das A, Iskratsch T and Krause S. *Biosensors and Bioelectronics*. Elsevier.

Mechanoregulation of vascular smooth muscle cells phenotype in atherosclerosis.

Swiatlowska P, Sit B, Feng Z, Marhuenda E, Zhou X, Xiao Q, Shanahan C, Yu C-H and Iskratsch T. *Journal of Molecular and Cellular Cardiology* vol. 173, 41-42. Elsevier.

Crosstalk between cell-matrix and cell-cell signaling in cardiomyocyte mechanosensing.

Marhuenda E, Iskratsch T and Xanthis I. *Journal of Molecular and Cellular Cardiology* vol. 173,. Elsevier.

The cardiomyocyte mechanical code Mechanical imprinting of talin interactions regulate cardiomyocyte mechanosensing and function.

Marhuenda E, Pandey P, Xanthis I, O'Neill C and Iskratsch T. *Journal of Molecular and Cellular Cardiology* vol. 173,. Elsevier.

Regulation of Cardiomyocyte Adhesion and Mechanosignalling Through Distinct Nanoscale Behaviour of Integrin Ligands Mimicking Healthy or Fibrotic ECM.

Iskratsch T, Hawkes W, Marhuenda E, Reynolds P, O'Neill C, PANDEY P, Wilson D, FREELEY M, HUANG D, Junquiang H, Gondarenko S, Hone J, Gadegaard N and Palma M. *Philosophical Transactions of The Royal Society B: Biological Sciences*. The Royal Society.

Pressure and stiffness sensing together regulate vascular smooth muscle cell phenotype switching.

Swiatlowska P, Sit B, Feng Z, Marhuenda E, Xanthis I, Zingaro S, Ward M, Zhou X, Xiao Q, Shanahan C, Jones GE, Yu C-H and Iskratsch T. *Science Advances* vol. 8, (15).

2021

Glioma stem cells invasive phenotype at optimal stiffness is driven by MGAT5 dependent mechanosensing.

Marhuenda E, Fabre C, Zhang C, Martin-Fernandez M, Iskratsch T, Saleh A, Bauchet L, Cambedouzou J, Hugnot JP, Duffau H, Dennis JW, Cornu D and Bakalara N. *Journal of Experimental and Clinical Cancer Research* vol. 40, (1).

Targeting Several Biologically Reported Targets of Glioblastoma Multiforme by Assaying 2D and 3D Cultured Cells.

Sixto-Lpez Y, Marhuenda E, García-Vazquez JB, Fragoso-Vazquez MJ, Rosales-Hernndez MC, Zacarías-Lara O, Méndez-Luna D, Gmez-Vidal JA, Cornu D, Norbert B and Correa-Basurto J. *Cellular and Molecular Neurobiology* vol. 42, (6) 1909-1920. *Springer Nature*.

2020

Matrix stiffness and blood pressure together regulate vascular smooth muscle cell phenotype switching.

Swiatlowska P, Sit B, Feng Z, Marhuenda E, Xanthis I, Zingaro S, Ward M, Zhou X, Xiao Q, Shanahan C, Jones GE, Yu C-H and Iskratsch T. *Cold Spring Harbor Laboratory*.

Glioma stem cells invasive phenotype at optimal stiffness is driven by MGAT5 dependent mechanosensing.

Marhuenda E, Fabre C, Zhang C, Martin-Fernandez M, Iskratsch T, Saleh A, Bauchet L, Cambedouzou J, Hugnot J-P, Duffau H, Dennis JW, Cornu D and Bakalara N. *Cold Spring Harbor Laboratory*.

2019

A novel 3D nanofibre scaffold conserves the plasticity of glioblastoma stem cell invasion by regulating galectin-3 and integrin-1 expression.

Saleh A, Marhuenda E, Fabre C, Hassani Z, Weille JD, Boukhaddaoui H, Guelfi S, Maldonado IL, Hugnot J-P, Duffau H, Bauchet L, Cornu D and Bakalara N. *Scientific Reports* vol. 9, (1). *Springer Nature*.